


At the Heart of **IMAGING**

cvi⁴²

Minimum Recommended Specifications

Standalone / Client



CPU

Intel Core i5



RAM

8 GB



VIDEO

Intel Integrated



DISK

60 MB/s
Effective
Transfer Rate



MAC OS


10.11 El Capitan



WINDOWS

Windows 7 SP1

Server



CPU

Intel Core i7
or Xeon



RAM

16 GB



DISK

60 MB/s
Effective
Transfer Rate



WINDOWS

Server 2008 R2

CT & 4D Flow

Higher Requirements



RAM

16 GB



VIDEO

Discrete nVidia or AMD
with at least 2GB vRAM



USERS

up to 3
Concurrent

*Please contact your
Circle advisor for
client-server
configurations
above 3
concurrent users*

Questions?

✉ support@circlecvi.com

☎ +1 (403) 338 1870